



A good friend to have around

Sami the Seal compressor nebulizer system

The Philips Respironics Sami the Seal pediatric compressor is shaped like a seal to appeal to children, their parents, and pediatricians, and to support aerosol therapy compliance. The compressor partners with our highly efficient SideStream nebulizers to provide fast and friendly treatment to pediatric patients in the home.

Character mask

For use with the Sami the Seal compressor nebulizer system, the Tucker the Turtle pediatric mask is also designed to appeal to children and support compliance. Tucker can be used with both the SideStream disposable and reusable nebulizers and has been proven to reduce aerosol deposition in and around the eyes.

PHILIPS
RESPIRONICS
sense and simplicity

Efficient, reliable and friendly

The Sami the Seal compressor nebulizer system combines power and reliability with a child-friendly design to create a nebulizer system that children will want to use. Sami is bottom weighted to help with stability and incorporates a docking station for the nebulizer between treatments. The system is available in two configurations – one for use with the SideStream disposable nebulizer only, and one with the SideStream disposable and SideStream reusable nebulizers.

Key advantages

- Reliable, economical and easy to use
- Bottom-weighted for stability
- Pediatric appeal
- Fast treatment times*
- Convenient handset dock



Ordering information

| Product description | Reorder number |
|---|----------------|
| Sami compressor with SideStream disposable nebulizer and Tucker mask | 1082456 |
| Sami compressor with SideStream disposable and reusable nebulizers and Tucker mask | 1082731 |
| Tucker the Turtle pediatric face mask (for use with SideStream reusable HS860 or SideStream disposable HS800) | 1044382 |
| SideStream disposable nebulizer | HS800 |
| SideStream reusable nebulizer | HS860 |
| Sami compressor replacement filter (5 pack) | 1083128 |



SideStream, Sami the Seal and Tucker the Turtle are trademarks of Respirationics, Inc. and its affiliates, all rights reserved.

© 2011 Koninklijke Philips Electronics N.V. All rights are reserved.

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication. CAUTION: US federal law restricts this device to sale by or on the order of a physician.

DD MCI 4104014 PN 1082451

SideStream disposable and reusable nebulizers

SideStream reusable and disposable nebulizers can be used with commonly prescribed respiratory medications and offer fast, efficient and economical drug delivery. Both nebulizers feature venturi technology that provides air flow in addition to that supplied by the compressor, resulting in faster drug delivery and shorter treatment times.*



Tucker the Turtle pediatric mask

The Tucker the Turtle child-friendly pediatric character mask is made of a soft, flexible material that contours to the face for comfort and fit. It also incorporates a vent design that directs medication away from the eyes, reducing aerosol deposition and minimizing the impact of long-term use.



* Data on file at Philips Respirationics

Product specifications

| | |
|-----------------------------------|---------------------------------|
| Electrical requirements | 120VAC-60 Hz |
| Power | 120 VA |
| Thermal overload protector | |
| Thermal fuse | Functioning temperature 130° C |
| Performance average flow rate | 6 LPM @ 10 psig |
| Maximum pressure | 29-43.5 psi |
| Size | 10.5 in L x 6.1 in W x 8.9 in H |
| Weight | 3.5 lbs. |
| Operating mode | Continuous |
| Warranty | 5 years |

Philips Healthcare is part of Royal Philips Electronics

Philips Respirationics International
Headquarters
+33 1 47 28 30 82
Philips Respirationics Asia Pacific
+65 6882 5282

Respirationics New Jersey, Inc.
5 Wood Hollow Road
Parsippany, NJ 07054
+1 800 345 6443
+1 724 387 4000
www.philips.com/respirationics